



Extend A-Constant conversion table

Optical / Immersion US
Contact US

	SRK/T A Constant	Hoffer Q pACD	Holloday 1 sf	a0	Haigis a1	a2
 BunnyLens	119.0	5.56	1.77	1.4	0.4	0.1
	118.4	5.23	1.44	1.03	0.4	0.1
 SeeLens	119.0	5.58	1.81	1.40	0.4	0.1
	118.5	5.25	1.49	1.05	0.4	0.1

1 IOL constant was evaluated using IOL master and the SRK/T formula, relations between constants

2 IOL constant was evaluated using US biometry and the SRK/T formula, relations between optical and us biometry

* It is recommended that surgeons personalize their IOL constant based on their surgical techniques and equipment, experience and post-operative results.