




Intensity A-Constant conversion table

Optical / Immersion US
Contact US

	SRK/T A Constant	Hoffer Q pACD	Holloday 1 sf	a0	Haigis a1	a2
 BunnyLens	118.4	5.15	1.39	0.928	0.4	0.1
	118.06	4.93	1.17	0.669	0.4	0.1
 SeeLens	118.4	5.15	1.39	0.928	0.4	0.1
	118.06	4.93	1.17	0.699	0.4	0.1
 Plate Haptic	117.45	4.65	0.82	0.346	0.4	0.1
	117.1	4.45	0.69	0.26	0.4	0.1

1 IOL constant was evaluated using IOL master and the SRK/T formula, relations between constants

2 IOL constant was evaluated using US biometry and the SRK/T formula, relations between optical and us biometry

* It is recommended that surgeons personalize their IOL constant based on their surgical techniques and equipment, experience and post-operative results.