

Hydrophilic Toric A-Constant conversion table

Optical / Immersion US
Contact US




Plate Haptic

117.7	4.86	1.02	0.448	0.4	0.1
117.3	4.61	0.77	0.184	0.4	0.1

1 IOL constant was evaluated using IOL master and the SRK/T formula, relations between constants

2 IOL constant was evaluated using US biometry and the SRK/T formula, relations between optical and us biometry

* It is recommended that surgeons personalize their IOL constant based on their surgical techniques and equipment, experience and post-operative results.