


# Hydrophilic Toric A-Constant conversion table

Optical / Immersion US  
Contact US

	Plate Haptic	117.7	4.86	1.02	0.448	0.4	0.1
		117.3	4.61	0.77	0.184	0.4	0.1

1 IOL constant was evaluated using IOL master and the SRK/T formula, relations between constants

2 IOL constant was evaluated using US biometry and the SRK/T formula, relations between optical and us biometry

\* It is recommended that surgeons personalize their IOL constant based on their surgical techniques and equipment, experience and post-operative results.